## Distortion correction of body diffusion weighted images using three point-Dixon method

S. Mori<sup>1</sup>, I. Muro<sup>1</sup>, H. Moriguchi<sup>1</sup>, T. Horie<sup>1</sup>, M. Honda<sup>1</sup>, T. Ogino<sup>2</sup>, M. Obara<sup>2</sup>, and Y. Imai<sup>1</sup>

<sup>1</sup>Radiology, Tokai university, Isehara, Kanagawa, Japan, <sup>2</sup>Philips Healthcare Asia Pacific, Shinagawa, Tokyo, Japan

Introduction: Significant distortion can often be observed in body diffusion weighted images (DWI) using single shot EPI. Distortion is caused by external magnetic field (B0) inhomogeneity. The degree of distortion is



Fig.1. Phantom (diameter 380mm)

proportional to B0. When a precise B0 map is obtained, distortion can be corrected using the B0 map. In this study, we demonstrate that a B0 map obtained using three point-Dixon method is useful to correct for distortion. Clinical images corrected using this method are also shown.

Method: A flow chart of our proposed method is shown in Fig.2. In this

image that produced the highest

as the final image. We performed

volunteer

phantom and

method, we applied five low pass filters to a B0 map. Only pass bands of these low pass filters were varied. Distortion correction was done using these five B0 maps. A final corrected image was selected from these five corrected images. We calculated correlation coefficients between each

abc

d e f

g h i




Fig.3. Sagittal MPR images ( reconstructed from axial images)

a,b,c: Images 170mm shifted from the center to the left

Images at the center

g,h,i: Images 170mm shifted from the center to the right

a,d.g: Images before distortion correction

b,e,h: Images corrected using B0 maps created using Sequence2

Images corrected using *c,f,l*: B0 maps created using Sequence3




Fig.2. A flow chart of distortion correction

scanning to test our proposed method. In phantom study, we scanned a column phantom filled with water (Fig.1) using eight channel phased array coils. The size of the phantom was 380mm x 380mm x 100mm. Details of scanning parameters were: FOV 400mm, slice thickness/gap 4/0mm, #slices 100, NEX 2~10, SENSE acceleration factor 2, image matrix size 256 x 256. We used three types of pulse sequences in this experiment. Parameters of these sequences are as follows. Sequence 1: single shot EPI TR/TE 7000/60msec; Sequence 2: three echo GRE (Dixon methods) #phase encoding (PE) steps 32, TR/TE 2500/3.2x3msec; Sequence 3: three echo- GRE (Dixon methods) #PE steps 128, TR/TE 2500/3.2x3msec. We created water and fat images, a B0 map and the combined water and fat image using the sequences 2 and 3. We performed distortion correction for 100 axial images. Two types of B0 maps were created using sequences 2 and 3. They differed in resolutions. Correction was done based on each B0 map. A sagittal multiplanar reformatted image was created using both of these corrected images. We made comparison on preciseness of the correction between these B0 maps. In in-vivo study, we scanned asymptomatic

volunteers. DWI of pelvis and spinal cords of the volunteers were scanned using eight channel phased array coils. b-value was set to 800s/mm<sup>2</sup>. We used the same pulse sequences as Sequences 2 and 3 described above in this experiment.

Results: In phantom study, there was no significant difference between the images corrected using the B0 map with 32 PEs and those corrected using the B0 map with 128 PEs. (Fig3) The correlation coefficients of both corrected images were increased by 0.01 or greater from those of the uncorrected images. In in-vivo study, significant improvements of distortion correction were

observed in both spinal cord and pelvis images.(Fig4) Discussion: In three point-Dixon method, good distortion correction is achieved using the sequence with 32 PEs. In this method, although prescan is required, it requires about a minute and a B0 map obtained from the prescan can be applied to 50 images. This method is quite useful in practice since a B0 map created using this method enables good distortion correction.

Acknowledgements: We thank Philips Electronics Japan for providing us 3 pt-DIXON(WIP).




Fig.4.In-vivo images

abcdefg

a: Pelvis DWI before distortion correction b: Combined water and fat images c: Pelvis DWI after distortion correction Proc. Intl. Soc. Mag. Reson Med 18 (2010)  $e: Sagittal\ MPR\ image\ without\ correction$   $f: Sagittal\ MPR\ image\ with\ correction$   $g: TSE\ conventional\ image\ f: Sagittal\ MPR\ image\ with\ correction$